


navigation
on the spot

World's first
*automatic
registration*



 **fiagon.tracey**



World's first
automatic registration with




 **fiagon.tracey**

www.fiagon.com/tracey

ENT navigation: safer, faster and simpler with Fiagon

- Fast and easy registration of the patient through registration via photograph with Fiagon.Tracey
- Ready-to-run – for inpatient and outpatient interventions
- Flexible instruments for challenging anatomical structures
- Precise navigation at the tip of the instrument
- Supports your workflow through maximum integration during the surgery
- For endoscopy and microscopy



Revolutionary navigation technology for routine and complex interventions

Wireless communication between the computer and measuring units and the iPad remote control



navigation on the spot.

Fiagon.Tracey offers you a completely new option for the central process of patient registration – registration via photograph.

The images are captured with the tablet and Fiagon.Tracey calculates the registration. As a result, the quality of the registration is independent of the user-specific instrument guidance – this is fast and easy for the surgeon and the patient. For this reason, it can be used during inpatient as well as outpatient procedures.



A worldwide innovation – Automatic registration

- Survey of the entire surface of the face through photo registration
- Fast and easy handling using the iPad remote control



Minimal footprint – 100% integration

- Hardware components reduced to a minimum
- Ready-to-run, for rapid inpatient and outpatient use



Flexible instruments navigated at the tip

- Electromagnetic navigation at the tip of the instrument
- Precise tracking of the flexible instruments in the surgical field