

Introducing

ENDURE™

ENDURE from Lumenis is a new revolution in laser fiber technology, enabling precise, smart tissue management across a wider range of procedures - including those with lower reimbursement. The ENDURE fiber perfectly balances dependability, quality, and flexibility with cost effectiveness. With this new innovation, surgeons can extend their capabilities to perform more surgeries with fewer outlays.


- Perform 5 surgeries per fiber
- High quality, dependable fiber with a tip that can be renewed during a procedure
- A clear aiming beam that helps you precisely target and treat only the tissue that you want
- Durable fiber with excellent energy transmittance
- Compatibility with flexible endoscopes for maximum access to the hard to reach areas

Compatible with all existing handpieces and accessories

Further Extending Your ENT Capabilities

Quality Beam Performance

Greater power with less beam divergence allows for maximal precision



ENDURE CO₂ Fiber

Part Number	AC-1148110
Specifications	2 meter long, 1.04mm outside diameter, 500µm inside diameter
Spot Size	295µm at fiber output
Number of Procedures	Performance tested for 5 typical ENT surgeries
Sterilization Method	PreVac Autoclave, Steam Sterilization
Accessories	For use with the ENDURE Sterilization Tray – (AC-1148700)

