

Introducing

# FiberLase™ STC



The FiberLase STC from Lumenis delivers exceptional dependability, quality, and flexibility.


A new ST connection system lets you benefit from FiberLase technology on your existing laser console. A two meter long, durable fiber gets the laser beam to exactly where you need it. With the FiberLase STC, you can attain better control of beam delivery with your laser for smarter tissue management.

- Durable fiber with excellent energy transmittance
- A tip that can be renewed during a procedure
- High quality, dependable fiber that's built to last
- A 2 meter long fiber
- Compatibility with flexible endoscopes for maximum access to the hard to reach areas

Further Extending Your ENT Capabilities

# Quality Beam Performance

Greater power with less beam divergence allows for maximal precision



	<b>FiberLase STC</b> (Part Number: AC-1059610)	<b>Competitor Fiber<sup>1</sup></b>
Fiber Length	2.0m	1.5m
Spot Size at Tip Output	295µm	320µm
Energy Transmission <i>Typical operating conditions</i>	>60%	Up to 50%
Outer Diameter	1.04mm	1.21mm
Renewable tip	YES	NO
Compatible Medical Grade Gas Types	Air, Helium	Helium

1. Data is based on published information in user manuals of other fiber manufacturers.

