

>>> PATIENT EYE PROTECTION



SPERIAN® EYESHIELDS

THE BEST TREATMENT REQUIRES THE BEST PROTECTION

GIVE THEM THE SAFE, EFFICIENT PROTECTION OF SPERIAN DISPOSABLE EYESHIELDS

For facial treatments, involving IPL, laser or microdermabrasion systems, Sperian's patented, hypo-allergenic design offers convenient, unobstructed access to nose bridge, temples and areas close to the eyes. For you and for your clients, Sperian disposable eyeshields are the safe, economical, easy-to-use protection solution.

Available for a variety of uses from simple physical barrier to full laser protection. **CE marked.**



SPERIAN® EYESHIELDS

All of Spierian's patented client eyeshields are easy to use, and allow easy access to nose bridge and crow's feet area, while providing protection against a variety of applications. All patient eyeshields are hypoallergenic, latex free and comfortable. Colorful packaging designs allow you to clearly distinguish the various models.



SPERIAN GPT IPL-AID™

- Provides protection from the wavelengths emitted from IPL devices hence blocking out the bright light commonly associated with IPL systems.

SPERIAN GPT LASER-AID®

- Provides protection from all lasers operating between 190-10,600nm with an optical density of 7. These shields also disperse the heat typically emitted from infrared lasers.

SPERIAN GPT DERM-AID®

- Provides a physical barrier for use during common dermatological procedures such as microdermabrasion, chemical peels and facials.

SPERIAN GPT LED-AID™

- Eyeshields are independently and internationally certified to withstand energy from LED devices of all wavelengths from 315-1400nm and low-powered (<1 watt) lasers. LED-Aid blocks out the bright light from red, yellow and blue LED's and IPL systems.

PART #	DESCRIPTION	SPECIFICATIONS	PACKAGING
31-7300	GPT Derm-Aid	For Non-Laser Use Only, CE Marked	50 Pair/Box
31-9300	GPT IPL-Aid	For Non-Laser Use Only, Opaque Inner Layer to Shield IPL, CE Marked	50 Pair/Box
31-8300	GPT Laser-Aid	Laser-Safe for Medical Lasers with Wavelengths Between 190nm – 11,000nm Protection Level: OD 7 @ 190-11,000nm EN 207: 315-1400 DIR L6, >1400-11000 D L4 + IR L3, GPT CE	24 Pair/Box
31-8400	LED-Aid	Safe for Use with LED's, Low-Powered Lasers (<1 Watt), IPL and Microdermabrasion Devices Protection Level: OD 4 @ 315-1400nm EN 207: >315-1400 DIR L4, GPT CE	50 Pair/Box

Spierian Eye & Face Protection, Inc., Laser Eyewear
 10 Thurber Blvd, Smithfield, RI 02917, USA | tel: 1.800.500.4739
 laserinfo@spierianprotection.com
 Christofsstrasse 5, D55116 Mainz, Germany | tel: +49.6131.9066497
 glendale-mainz@spierianprotection.com

SEF.225 9/09

ALSO AVAILABLE

SPERIAN IPL EYEWEAR

Spierian offers a wide variety of frame styles for IPL applications. Available with light, medium or dark shades to protect you from the radiation and bright lights associated with IPL devices.

GPT XC®: Features Multi Material Technology (MMT™) that cushions every area where the spectacle comes in contact with your face. The result is 180 degrees of comfort. Optional Rx insert is available.



NEW! SPERIAN MILAN™: Fashion-inspired, features form-fitting lenses, and multiple adjustments in temple and nose-bridge areas that allow you to personalize the fit to your face.



GPT LGF™: This is a wrap-around spectacle that also fits over many prescription eyewear styles.



ORDERING INFORMATION

SHADE	ITEM #	FRAME & SHADE TYPE
Light	31-80133	GPT XC, Light Shade
	31-60133	Spierian Milan, Light Shade
	31-9623	GPT LGF, Light Shade
Medium	31-80141	GPT XC, Medium Shade
	31-30141	GPT LGF, Medium Shade
Dark	31-80134	GPT XC, Dark Shade
	31-9624	GPT LGF, Dark Shade

Distributed by