



Legend Pro+

Dare to imagine



4 in 1 platform

Personalized aesthetic solutions that meet your patients' expectations and more

Imagine a treatment with immediate results,
long-term efficacy, minimal pain and downtime,
and no need for anesthesia...

Sounds too good to be true?

Legend Pro+™ is all that, and more.

The platform with the power of synergized technologies

Legend Pro+™ offers endless options for personalized head-to-toe treatment plans. From muscle activation, to skin remodeling and resurfacing to address the specific needs and desires of each patient.

With **Legend M²™** combination face treatment and Legend combination of muscle toning and skin tightening treatment, you can get optimal results and maximize patient satisfaction.

- 4 technologies in 1 platform
- Treats the skin from inside out
- Minimal pain and downtime
- Immediately visible results
- Proven combined therapy protocol



Make your patients' wishes come true while ensuring their safety and long-term satisfaction

Each patient is a world of their own; complete, unique and singular.

Legend Pro+™ allows them to enjoy that individuality by enhancing their best characters.

The Legend Pro+™ platform offers a complete tool kit for tailoring the right treatment to each of your patients, with the flexibility of treating a wide range of skin types, and imperfections.

Treatment with **Legend M²™** reaches multiple layers of the skin, introducing a wide range of solutions for skin resurfacing and inducing regeneration, as well as wrinkle reduction. Combination of muscle toning and TriPollar[®] RF provides solutions from lifting tightening and treating the skin from within approach.

Testimonials

"The TriPollar[®] RF and the Hybrid Energy allow me to customize the treatments for my patients and to get optimal results. My patients are happy going through almost painless treatment, and the treatment time is shortened without the need for numbing. These unique technologies allow me to treat multiple indications and get exceptional results."

Matteo Tretti Clementoni, MD – Laserplast Sagl, Milan, Italy

"Hybrid Energy Technology treatments result in a relatively major dermal remodeling effect while minimizing visible epidermal effect. This treatment modality is very suitable for delicate and sensitive skin types which are prone to pigment changes"

Michael H. Gold, M.D., Nashville, Tennessee, USA

"With combination of RF and dynamic muscle activation, the efficacy of the treatment of cellulite appearance, skin tightening, and body shaping is supposed to increased when comparing to the treatment with RF alone."

Wanitphakdeedecha R. Treatment of abdominal cellulite and circumference reduction with radiofrequency and dynamic muscle activation. Journal of Cosmetic and Laser Therapy, April 2015

Seeing is Beliveing



TRIPOLLAR



DMA

Lumenis Collection



Before

After 5 treatments

Courtesy of Alex Levenberg, MD



HYBRID ENERGY



Before

After 3 treatments

Courtesy of Matteo Tretti Clementoni, MD



HYBRID ENERGY



Before

After 4 treatments


By encouraging and enhancing the natural processes of the body, Legend Pro+™ unveils the best in each patient.

From resurfacing of the upper layers of the skin, promoting the regeneration of collagen and elastin fibers in the deeper dermis to toning and firming the SMAS layer, firming and toning of muscles in the treatment area, Legend Pro+™ offers a complete toolkit of non-invasive and minimally invasive technologies for creating tailor-made treatments and personalized results for your patients.



TriPollar®


Triggering regeneration of collagen and elastin fibers, by delivering focused RF energy in a safe and effective manner, inducing controlled heating of the deep dermis.



DMA™

Dynamic Muscle Activation

In motion applicator activating an electrical current for non-invasive SMAS stimulation for **toning and lifting of the SMAS** layer to create a "Face lift" effect and toning and shaping of muscle groups in the treatment area.




Hybrid Energy™

Activating the body's natural healing processes by applying RF assisted thermal ultra-thin electrodes ablation, followed by galvanic current to the tissue to create a non thermal effect creating micro channels to induce controlled micro wounds, thus **encouraging the generation of new and healthy skin**.


Smart Tip Design

Unique smart tip design enabling even distribution of RF energy that flows between the ultra-thin electrode and the disposable tip's surface, thus **reducing the risk of arching and undesired side effects**.



TriFractional™

Applying RF electrode ablation with controlled volumetric heating to achieve **skin resurfacing and regeneration** via wound healing process, resulting in skin texture improvement and refinement.



Built for perfection for both patients and physicians

Legend Pro+™ includes unique design features that are created to ensure the safety and comfort of the treatment to patients and physicians.

RF flow between ultra-thin electrodes and tip's surface, **ensuring safer RF administration**

Real-time in-motion milliseconds precision of temperature reading, to increase confidence in clinical results by **ensuring optimal-temperature throughout the treatment**

Several disposable tips configurations, differing in length and density, **to ensure flexibility for creating the desired clinical effect**

Two technologies working simultaneously in a single applicator reaching both **skin and muscle** for improved outcome

Patient biofeedback button allowing patients to **communicate their level of comfort during the treatment**

Ergonomically designed to allow a **close fit to the patient's curves**, as well as to ensure **a comfortable treatment by the physician**



The proprietary technologies in **Legend Pro+™** are scientifically proven through clinical research and market experience.

Legend M2™ clinically proven and tested combination therapy

During the **Legend M2™** treatment the fibroblasts undergo two different triggers, minimally invasive, and non-invasive. The dual effect creates a synergetic result.

Legend M2™ unique synergized treatment combines the **macro** volumetric heating of tissue with TriPollar® RF and the precise **micro** heating of focal points with HE-A RF- assisted ultra-thin electrodes to turn your patients' wishes into reality.

Courtesy of Mark Rubin, M.D



Before

After 3 treatments

Legend M2™ treatment validated through clinical experience, as published in the article by Steven D. Shapiro "VoluDerm™* Micro-Needle Technology for Treating Skin Laxity and Wrinkles-Initial Clinical Experience".

Jacobs Journal of Experimental Dermatology, January 2015

* Hybrid Energy Program A



TRIPOLLAR HYBRID ENERGY

Muscle stimulation and TriPollar® combination

Legend Pro+™ combination treatment of muscle toning and skin tightening offers a new approach of treating the skin from within.

Volumetric heating of the skin with **TriPollar®** RF, leads to collagen remodeling and skin tightening followed by **DMA™ muscle stimulation**, leading to lifting and toning of the SMAS layer.

For treatment areas other than the face, the **TriPollar® RF volumetric heating together with DMA™ muscle activation**, lead to skin tightening and toning and shaping of muscle groups in the treatment area.

Lumenis collection



Before

After 6 treatments

Before

After 6 treatments



TRIPOLLAR

DMA

Your journey with Lumenis has just begun

We want your business to grow with Lumenis. That's why we have put together a support package that begins with the first conversation and will stay with you as your business begins to see the benefits. This ensures you have all the training you need and receive the most up to date clinical materials. Our direct to consumer marketing activities support the brand, keeping Legend Pro+™ the treatment of choice when it comes to wrinkle, resurfacing, skin tightening and muscle toning treatments.



Lumenis Practice Locator

Thousands of patients search for treatment providers every month. Help them find your clinic by registering your practice on the Lumenis Practice Locator. We run campaigns through influencers, social media and viral videos to promote this directory.



Lumenis online

Consumers have many questions about aesthetic. Our patient-facing websites will give them the re-assurance and confidence they need to move forward with the treatments.



Team Training

We install the system and train your staff, so that you are able to assess skin types and conditions, creating the perfect treatment plan for each individual patient. Our expert clinical training means they will handle consultations and treatments with confidence.



Your Most Valuable Resource

Your dedicated representative is an email, text, or phone call away. We want you to feel part of the Lumenis family, so please feel free to contact us any time. We'll always do our best to get back to you within 24 hours.



Brand Your Clinic

As a Lumenis client you will have access to the Lumenis Physician Portal. Here, you can find materials such as brochures and posters. You will also have access to a wealth of clinical and product information.



Extensive Service Network

We are committed to providing high quality, comprehensive service for medical equipment repair and maintenance. Lumenis service provides you with peace of mind, so you can focus on patient care, rather than on the equipment.

Specifications

Legend Pro+™

Parameter	Value/Data
Input voltage	100-240V~, 50/60 Hz, max 2.2A
Mode of operation	Hybrid Energy™, TriFractional™, TriPollar®, DMA™
Preset Parameters	Yes
Maximum output power RF	50 Watts @ 200 Ohm
Maximum output voltage	100±10 Vrms
Output frequency RF	1MHz
Platform Specifications	
Weight	~30 Kgs / 66 lbs
Dimensions	L45 x W45 x H110 cm

TriPollar® and DMA™ Applicators

TriPollar®	Number of electrodes	Temperature sensor
Applicator 1	6	Integral
Applicator 2	3	Integral
Applicator 3	3	External

Hybrid Energy™ Applicator

Hybrid Energy™ tips	Ultra-thin electrodes array	Length	Pulse count
Gen36	6x6	0.6 mm	Up to 800
Gen12	6x2	0.6 mm	Up to 800
Gen100	10x10	0.6 mm	Up to 800
Gen36L	6x6	1.0 mm	Up to 800
TriFractional™ tips	Pins array	Electrode length	Pulse count
H7X7	7x7	0.2 mm	Up to 800

Warnings and risks:

Legend Pro+™ is contraindicated for patients with pacemaker, defibrillator, or any implanted electronic device, metal implants in the treatment area. Improper use of the System could result in possible side effects.

Although these effects are rare and expected to be temporary, side effects may appear either at the time of treatment or shortly after. Side effects may include any of the following: prolonged or significant pain damage to natural skin texture (blister, burn), excessive skin redness (erythema), excessive swelling (edema), fragile skin bruising, excessive itching, change of pigmentation (hyper-pigmentation or hypo-pigmentation), scarring transient skin break-out such as acne and pimples. For complete list of contraindications, please refer to the User Manual.

About Lumenis®

Lumenis is a global leader in the field of minimally invasive clinical solutions for the surgical, ophthalmology and aesthetic markets, and is a world-renowned expert in developing and selling energy-based technologies, including laser, intense pulsed light (IPL) and radio frequency (RF).

For over 50 years, Lumenis' ground-breaking products have redefined medical treatments and have set numerous technological and clinical gold standards.

© 2018 All Rights Reserved. The Lumenis Group of Companies. Lumenis, its Logo, Legend Pro+, Hybrid Energy, DMA, TriPollar and TriFractional are trademarks or registered trademarks of The Lumenis Group of Companies. Specifications are subject to change without notice. Hybrid Energy, DMA, TriPollar and TriFractional are proprietary technologies of The Lumenis Group of Companies.

 **Pollogen Ltd.,**
6 Kaufman St.,
P.O.Box #50320
Tel Aviv 6801298, Israel

Lumenis LTD
Yokneam Industrial Park
Hakidma Street 6,
P.O.B. # 240
Yokneam 2069204, ISRAEL
T +972 4 959 9000

EMEA
Dreieich Dreieichenhain, Germany
T +49 6103 8335 0
F +49 6103 8335 300

Roma, Italy
T +39 06 90 75 230
F +39 06 90 75 269

Hertfordshire, UK
T +44 20 8736 4110
F +44 20 8736 4119

JAPAN
Tokyo, Japan

T +81 3 6743 8300
F +81 3 6743 8301

ASIA PACIFIC
Beijing, China
T +86 10 5737 6677
F +86 10 5737 6767

Gurgaon, India
T +91 124 421 07 95

Kowloon, Hong Kong
T +852 217 428 00
F +852 272 251 51

© Lumenis Ltd 2018