

# ResurFX™

True Fractional Non-Ablative Skin Resurfacing



# Advanced. Intuitive. Affordable.



ResurFX™ is an advanced 1565nm fiber laser system for True fractional non-ablative skin resurfacing. The ResurFX sophisticated CoolScan scanner, enables you to offer better outcomes for a variety of patients, with NO disposables.

Clinical Results	Before	After
Courtesy of Matteo Tretti Clementoni, MD.		
Courtesy of Mitchel P. Goldman, MD.		
Courtesy of Matteo Tretti Clementoni, MD.		















*“As a dermatologist, I’m always looking for laser treatments that strike the right balance between efficacy and patient comfort. With ResurFX, we have a solution that improves facial and body skin tone and texture with minimal discomfort and downtime.”*

—  
Dr. Arielle N. B. Kauvar, MD,  
Director New York Laser & Skin Care

# CoolScan™ for Better Patient Outcomes

## The only true fractional non-ablative technology

- State-of-the-art CoolScan scanner for customizable treatment
- Patent pending algorithm ensures precise, homogeneous and safe energy delivery
- Efficacious and fast, with up to 20% coverage in one pass
- The only one with continuous contact cooling, for safe treatment and patient comfort
- Cost effective with NO disposables




# Specifications

## ResurFX™

Wavelength	1565 nm
Pulse Energy	10-70 mJ per micro-beam
Pulse Characteristic	CoolScan™ Scanner
Pattern Shapes	Line, square, rectangle, circle, donut, hexagon
Beam Density	Up to 500 micro-beams/cm <sup>2</sup>
Repetition Rate	0.5-2 Hz
Tips	Sapphire cooled tip & Precision tip
Spot Size	Up to 18 mm
Cooling	Continuous contact cooling
Dimensions	51 x 60 x 42 cm (20.1 x 28.6 x 16.55 inches)
Weight	22 kg / 48.5 lbs
Electrical requirements	100 – 240 VAC ±10%, 6A max., 50/60 Hz, single phase, dedicated line

© 2014 All Rights Reserved. The Lumenis Group of Companies. PB-1007789 Rev B. ResurFX and CoolScan are trademarks or registered trademarks of the Lumenis Group of Companies. Specifications are subject to change without notice. Lumenis is ISO 13485:2003 certified. ResurFX is CE pending. | ResurFX is pending FDA clearance.



For more information:  
[www.lumenis.com](http://www.lumenis.com)

### AMERICAS

San Jose, CA, USA  
T +1 408 764 3000  
F +1 408 764 3999

### EMEA

Dreieich Dreieichenhain,  
Germany  
T +49 6103 8335 0  
F +49 6103 8335 300  
Roma, Italy  
T +39 06 90 75 230  
F +39 06 90 75 269  
Hertfordshire, UK  
T +44 20 8736 4110  
F +44 20 8736 4119

### JAPAN

Tokyo, Japan  
T +81 3 6743 8300  
F +81 3 6743 8301

### ASIA PACIFIC

Beijing, China  
T + 86 10 5737 6677  
F +86 10 5737 6767  
Gurgaon, India  
T + 91 124 422 07 95  
Kowloon, Hong Kong  
T + 852 217 428 00  
F + 852 272 251 51